

OtoClear® Adapter Wand Instructions

For Use with Portable Waterpik® (Reorder #7215)



1 Push the adapter wand into the positioning wheel.



2 Twist the OtoClear® tip on the luer lock fitting of the adapter wand. (Be sure to tighten the OtoClear® tip securely).

NOTE: The wand can be cleaned as needed depending on usage rate. Use a soft cloth and mild, non-abrasive cleaner.



Contents: 3 Adapter Wands.



R_x ONLY



©2019 Bionix Development Corp.
5154 Enterprise Blvd., Toledo, Ohio 43612
Phone: 800.551.7096 | Fax: 800.455.5678
Email: bionix@bionix.com | bionix.com

*OtoClear® Tips - U.S. Patent #6,706,023.
Waterpik® is a registered trademark of Water Pik, Inc.
The OtoClear Safe Irrigation Tip is not affiliated with
Waterpik Technologies.

RM95-0059 | Rev. B